



### **Specifications**

Dimensions:	15 ½ W x 18 L x 14 H
Weight:	40 lbs. (without battery) 42 lbs. (with battery).
Electrical:	50/60 Hz 110-240 Vac, 24 volt DC power port.
Laser Classification:	Class 1 Laser Product. Compliance with FDA HHS 21 CFR 1040.10 and IEC 60825-1.
Imaging Media:	Phosphor Storage Plates, all brands. All standard and custom sizes.
System Includes:	Power Supply, DC Power Adapter and USB Cable.
Software:	Compatible, third party software recommended (not included).
Part Numbers:	D5000K-Q (ScanX Quantum w/o battery) D5000K-QB (ScanX Quantum w/ battery)

### **Computer Recommendations**

CPU Speed:	3.0 GHz Intel CORE2
Operating System:	Microsoft Windows XP Professional (for a 32-bit Intel processor with Service Pack 2 or later or for a 64-bit Intel x64 processor) Microsoft Vista Business, Enterprise, or Ultimate (for a 32-bit Intel processor with Service Pack 1 or later or for a 64-bit Intel x64 processor with Service Pack 1 or later).
RAM:	2 GB
CD ROM Drive:	CD/RW/DVD
Hard Drive:	500 GB
Optical Drive:	Optical drive that is CD Rom compatible
Monitor:	SVGA 20", 1024x768 or higher resolution, contrast ratio 450:1, .22 dot pitch recommended.
Video Display Adapter:	32 MB RAM
Keyboard / Mouse:	Standard keyboard and mouse
USB Port/Version:	USB 2.0 or later