

NEW

VET-TEC 2000

WALK- ON SCALE

WITH ACCU-LOCK II



TECHNIDYNE CORPORATION

www.technivet.com

1-800-654-8073

Vet- Tec 2000

The Perfect Scale for Your Practice

New ACCU-LOCK II System assures accuracy with moving animals

The design of the Vet-Tec 2000 is based on a newly developed single chip computer. This powerful chip is programmed in "C++" by a former Bell Labs programmer. It will display a more accurate weight on a moving animal than other scales because it computes the animal's true weight based on 300 readings taken over a 10 second period. The operation is completely automatic. Just walk the animal on the scale and the weighing process begins-it is not necessary to push any buttons.

Extra Bright Display, Easy to Read

New generation LED display offers exceptional brightness at wide viewing angles. Ideal for use in all lighting conditions.

User Friendly

Once turned on, it is always ready to weigh the next animal. Just walk the animal on the Low-Profile Platform (2"H) and the weighing process starts automatically, immediately displaying the weight without the need to push any buttons. If you forget the weight after the animal gets off - just press the LAST READ button and the weight will reappear.

Self - Leveling Platform

Since most floors are not perfectly level, the Vet-Tec 2000 was designed to automatically compensate for a non-level floor which prevents rocking.

2 Week Free Trial

Satisfaction guaranteed.

If during the first two weeks you are not satisfied, return the scale for a full refund.

3 Year Warranty

Unlike other scale companies, Technidyne offers a full 3 year warranty which includes parts and labor. Technidyne is the manufacturer of the scale and therefore will always be able to service the scale - no matter how old it is! Other Veterinary suppliers contract out to different scale manufacturers, which change from time to time - depending on who is the lowest bidder. They can then no longer service their older scales.

Special Features Include:

- Low Profile - 2"H
- Large Platform - 42"w x 22"d
- Heavy Gauge Stainless Steel Construction
- Rubber Mat Included
- Non-Slip / Self-Leveling Feet
- No Batteries Required
- Fully Assembled
- Converts Lbs/Kgs

Specifications:

- Accuracy- 0.1% at 300 lbs.
- Resolution- Increments of 0.1lb
- Max Load wt- 450 lbs
- Power Consumption - 5 watts
- Shipping wt- 45 lbs
- Measures up to 350 lbs

Technidyne Corporation

2190 Route 9, Unit D2,
Toms River, NJ 08753

www.technivet.com

1-800-654-8073